

# SOFT FORM® MATERNITY SUPPORT BELT

Model No: 31-570

Sizes: Small - Large  
Color: Beige

US Patent D437,416

Size according to pre-pregnancy dress size.	
Size	Fits
Small	2 - 6
Univ (Med)	6 - 14
Large	16 - 20

## MATERIAL CONTENT:

40% polyester, 25% nylon, 15% polyurethane, 10% neoprene, 5% lycra, 5% cotton. Latex free.

800-327-4110

## INDICATIONS:

Helps relieve the minor lumbar and abdominal pain and discomfort associated with pregnancy.

- Muscle Strain
- Lower Back Pain and Discomfort
- Abdominal Support

## BENEFITS:

Comfortably supports the lumbar and abdominal regions to help provide relief of lower back pain and discomfort. Unique multi-layered soft laminate provides excellent support and comfort. Ventilated, latex-free elastic compression side panels provide compressive support. Lightweight flexible stays help keep the support in place and prevent the belt from rolling. Can be worn inconspicuously under most clothing.





## PRODUCT FEATURES:

- **Abdominal Lift Attachment** provides abdominal lift with no undesirable pressure. Expands to accommodate growth during various stages of pregnancy.
- **Ventilated Elastic Back and Side Panels.** Latex-free ventilated elastic provides comfortable and cool support compression.
- **Unique Laminate Material** combines three materials: durable outer in an attractive beige color, support foam base and soft cotton inner lining.
- **Easy Hook and Loop Closure** allows for easy application.
- **Two Anti-Roll Stays** are precisely placed to help keep the support in place and prevent rolling.
- **Latex Free**

## HOW TO APPLY:

1. Wrap the Abdominal Lift Attachment around your waist. Position the pad so it is centered and under the abdomen. Wrap the elastic strap around your body and attach it to the pad. Both sides of the strap can be adjusted for a good fit. Ensure the abdominal pad is properly positioned but does not interfere with walking or sitting. The fit of this pad can expand as your baby grows to ensure proper fit. Simply fasten farther away from the center of the pad on the Abdominal Lift Attachment.
2. Attach the belt to the pad. With the belt centered on the lower back, pull both elastic sides until the desired tension and attach the closure to the abdominal pad. Ensure the "hook" ends of the belt have completely engaged to the pad to prevent skin irritation or snagging of clothing.
3. To help keep the support in place and keep from sliding down, stretch the top strap over your stomach and attach to the opposite side.
4. For even greater lift and support, with the support on, lay on your back with your knees up. Unfasten the belt and retighten the elastic strap and reattach closures.



## CARE INSTRUCTIONS:

Fasten all closures. Hand wash in cool water with mild soap. Remove excess water and air dry away from heat. Wash as needed.