

# CARPALMATE® WRIST SUPPORT

**Model No: 22-140**

**Size: Universal**

**Color: Black or Beige**

US Patents D463,566; 6,517,501;  
& D468,434

**800-327-4110**

## USE:

Helps prevent and provide symptomatic relief of Carpal Tunnel Syndrome. Provides support when performing tasks such as typing, data entry, light assembly, or any repetitive tasks. Comfortable for all day wear.



The unique design and angle of the CarpalMate® elevate the wrist and alleviate pressure on the Median nerve, which is the cause of pain and discomfort in Carpal Tunnel Syndrome (CTS). The snug fit of CarpalMate® at the heel of the palm limits hand and wrist motions, such as drooping and pinching, which aggravate CTS.

## COMFORT BENEFITS:

- **Compression-formed support base** is molded to conform to the wrist and palm area. The CarpalMate® provides support to the wrist by maintaining it in a neutral position. The Palmar pad provides a resting point for the hand and inhibits thumb and little finger opposition.
- **Low Profile Support** provides for full finger and hand dexterity while controlling thumb and little finger opposition – no uncomfortable thumb loops or bulky straps.
- **Rigid stay** maintains the wrist in orthopedically correct angle, alleviating pressure on the Median Nerve.
- **Wide elastic strap** provides adjustable compressive support over the wrist area.

## WEARING INSTRUCTIONS:

Unfasten closure. Open strap enough to pass hand through elastic strap and brace. Lay your hand in the brace with the heel of your hand against the front of the brace. Form a fist, allowing your fingertips to further position the brace against the palm of your hand. Use other hand to comfortably tighten the elastic strap..

## CARE INSTRUCTIONS:

Hand wash in cold water with mild detergent. Pat dry with towel to remove excess water. Air dry away from heat.



Website: **[www.flaorthopedics.com](http://www.flaorthopedics.com)**

©2004 FLA ORTHOPEDICS, INC.